

Calibration of and Preliminary Results for a NUMI Beam Coincident Neutron Detector in the MINOS Tunnel

Zachary Pierpiont

Penn State University

June 30, 2008



Mentor: Andrew Sonnenschein

Overview

- Introduction
 - Coupp
 - My Detector
- Calibration
 - Energy Calibration
 - Separation Algorithm
- Data Acquisition
- Results
- Future Work

COUPP

- Who? Fermilab, University of Chicago, Indiana University South Bend
- What? A bubble chamber dark matter detector
- Where? Currently in the MINOS tunnel, but in the future, SOUDAN
- Why? Excellent intrinsic gamma rejection
Good sensitivity to both SD and SI interactions
Low cost




My Detector

The detector is a 2 in. x 2in. liquid scintillator hooked up to a PMT.


Data was read via a digital oscilloscope.

The goal of this detector was to study background radiation in the MINOS tunnel, specifically beam coincident neutrons.




Gamma Energy Calibration

Histogram of energy from Cs-137 source in arbitrary units



Neutron Energy Calibration




Arb. Units	Neutron Energy MeV
1	1.7
5	5.0
10	8.2
50	30
100	58

$$N = .56 * A + 1.9$$

ICARUS Collaboration, "Calibration of BC501A Liquid Scintillator cells with monochromatic neutron beams." 13 March 1998

Gamma/Neutron Separation




Algorithm: Tail to Total Ratio

Q_{tot} = Integral of 130 bins from beginning of peak


Q_{del} = Integral of 130 bins starting 30 bins after beginning of peak

The variable used for separation is Q_{del}/Q_{tot}


Gamma/Neutron Separation



Tunnel Setup




Data



Coincident Neutron Energies

Histogram of Probable Neutron Energies



27 Probable Neutrons
375,000 Beam Pulses
Highest Energy: 19.32 MeV

Future Work

Take more tunnel data:

- Attempt to develop a more accurate coincident neutron spectrum

Improve energy calibration:

- Determine if Cs-137 peak was full absorption or Compton scattering peak

- Explore other neutron calibration methods

Better separation algorithm:

- Pulse Shape Analysis, a three term exponential fit

Shielding studies

- Use this detector in conjunction with different shields to test their effectiveness

Acknowledgments and Thanks

Andrew Sonnenschein

Erik Ramberg

ICARUS Collaboration

Nuclear Innovation Unit at CIEMAT

Penn State University

Department of Energy

Any Questions?